



IFT74 Antibody, FITC conjugated

Product Code	CSB-PA846657LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96LB3
Immunogen	Recombinant Human Intraflagellar transport protein 74 homolog protein (375-600AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Intraflagellar transport protein 74 homolog (Capillary morphogenesis gene 1 protein) (CMG-1) (Coiled-coil domain-containing protein 2), IFT74, CCDC2 CMG1
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	IFT74
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.