



# PARP8 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA847623LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q8N3A8
<b>Immunogen</b>	Recombinant Human Poly [ADP-ribose] polymerase 8 protein (628-854AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Poly [ADP-ribose] polymerase 8 (PARP-8) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 16) (ARTD16), PARP8
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Epigenetics and Nuclear Signaling
<b>Target Names</b>	PARP8