



PARP8 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA847623LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q8N3A8 |
| Immunogen | Recombinant Human Poly [ADP-ribose] polymerase 8 protein (628-854AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Poly [ADP-ribose] polymerase 8 (PARP-8) (EC 2.4.2.30) (ADP-ribosyltransferase diphtheria toxin-like 16) (ARTD16), PARP8 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | PARP8 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |