



ompC Antibody, HRP conjugated

Product Code	CSB-PA848941LB01EOD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8XE41
Immunogen	Recombinant Escherichia coli O157:H7 Outer membrane protein C protein (22-367AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli O157:H7
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein C (Outer membrane protein 1B) (Porin OmpC), ompC
Immunogen Species	Escherichia coli O157:H7
Research Area	Others
Target Names	ompC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.