



MICAL1 Antibody, HRP conjugated

Product Code	CSB-PA851555LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TDZ2
Immunogen	Recombinant Human [F-actin]-methionine sulfoxide oxidase MICAL1 protein (806-1067AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	[F-actin]-monooxygenase MICAL1 (EC 1.14.13.225) (Molecule interacting with CasL protein 1) (MICAL-1) (NEDD9-interacting protein with calponin homology and LIM domains), MICAL1, MICAL NICAL
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MICAL1