



# MICAL1 Antibody, Biotin conjugated

|                            |                                                                                                                                                                                              |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA851555LD01HU                                                                                                                                                                           |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                |
| <b>Uniprot No.</b>         | Q8TDZ2                                                                                                                                                                                       |
| <b>Immunogen</b>           | Recombinant Human [F-actin]-methionine sulfoxide oxidase MICAL1 protein (806-1067AA)                                                                                                         |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                       |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                        |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                        |
| <b>Form</b>                | Liquid                                                                                                                                                                                       |
| <b>Conjugate</b>           | Biotin                                                                                                                                                                                       |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                                             |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                                     |
| <b>Isotype</b>             | IgG                                                                                                                                                                                          |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                   |
| <b>Alias</b>               | [F-actin]-monooxygenase MICAL1 (EC 1.14.13.225) (Molecule interacting with CasL protein 1) (MICAL-1) (NEDD9-interacting protein with calponin homology and LIM domains), MICAL1, MICAL NICAL |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                         |
| <b>Research Area</b>       | Tags & Cell Markers                                                                                                                                                                          |
| <b>Target Names</b>        | MICAL1                                                                                                                                                                                       |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                  |