



MKL1 Antibody, Biotin conjugated

Product Code	CSB-PA853375LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q969V6
Immunogen	Recombinant Human MKL/myocardin-like protein 1 protein (674-812AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MKL/myocardin-like protein 1 (Megakaryoblastic leukemia 1 protein) (Megakaryocytic acute leukemia protein) (Myocardin-related transcription factor A) (MRTF-A), MKL1, KIAA1438 MAL
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MKL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.