



FCRL5 Antibody, HRP conjugated

Product Code	CSB-PA853494LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96RD9
Immunogen	Recombinant Human Fc receptor-like protein 5 protein (419-553AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fc receptor-like protein 5 (FcR-like protein 5) (FcRL5) (BXMAS1) (Fc receptor homolog 5) (FcRH5) (Immune receptor translocation-associated protein 2) (CD antigen CD307e), FCRL5, FCRH5 IRTA2
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	FCRL5