



# FCRL5 Antibody, FITC conjugated

|                            |   |
|----------------------------|---|
| <b>Product Code</b>        | CSB-PA853494LC01HU  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.   |
| <b>Uniprot No.</b>         | Q96RD9  |
| <b>Immunogen</b>           | Recombinant Human Fc receptor-like protein 5 protein (419-553AA)  |
| <b>Raised In</b>           | Rabbit  |
| <b>Species Reactivity</b>  | Human   |
| <b>Form</b>                | Liquid  |
| <b>Conjugate</b>           | FITC  |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4  |
| <b>Purification Method</b> | >95%, Protein G purified  |
| <b>Isotype</b>             | IgG   |
| <b>Clonality</b>           | Polyclonal  |
| <b>Alias</b>               | Fc receptor-like protein 5 (FcR-like protein 5) (FcRL5) (BXMAS1) (Fc receptor homolog 5) (FcRH5) (Immune receptor translocation-associated protein 2) (CD antigen CD307e), FCRL5, FCRH5 IRTA2 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)  |
| <b>Research Area</b>       | Immunology  |
| <b>Target Names</b>        | FCRL5   |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.   |