



CSF3R Antibody, Biotin conjugated

Product Code	CSB-PA860321LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99062
Immunogen	Recombinant Human Granulocyte colony-stimulating factor receptor protein (117-337AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granulocyte colony-stimulating factor receptor (G-CSF receptor) (G-CSF-R) (CD antigen CD114), CSF3R, GCSFR
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CSF3R
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.