



KCNN3 Antibody, Biotin conjugated

Product Code	CSB-PA880079LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UGI6
Immunogen	Recombinant Human Small conductance calcium-activated potassium channel protein 3 protein (240-426AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Small conductance calcium-activated potassium channel protein 3 (SK3) (SKCa 3) (SKCa3) (KCa2.3), KCNN3, K3
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	KCNN3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.