



MSRA Antibody, FITC conjugated

Product Code	CSB-PA883456LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UJ68
Immunogen	Recombinant Human Mitochondrial peptide methionine sulfoxide reductase protein (24-235AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mitochondrial peptide methionine sulfoxide reductase (EC 1.8.4.11) (Peptide-methionine (S)-S-oxide reductase) (Peptide Met(O) reductase) (Protein-methionine-S-oxide reductase) (PMSR), MSRA
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MSRA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.