




GIMAP4 Antibody

| | |
|----------------------------|--|
| Product Code | CSB-PA885732LA01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9NUV9 |
| Immunogen | Recombinant Human GTPase IMAP family member 4 protein (5-101AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA, IF; Recommended dilution: IF:1:50-1:200 |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | GTPase IMAP family member 4 (Immunity-associated nucleotide 1 protein) (IAN-1) (hIAN1) (Immunity-associated protein 4), GIMAP4, IAN1 IMAP4 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Immunology |
| Target Names | GIMAP4 |

Image



Immunofluorescence staining of HepG2 cells with CSB-PA885732LA01HU at 1:100, counter-stained with DAPI. The cells were fixed in 4% formaldehyde, permeabilized using 0.2% Triton X-100 and blocked in 10% normal Goat Serum. The cells were then incubated with the antibody overnight at 4°C. The secondary antibody was Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.