



# GIMAP4 Antibody, Biotin conjugated

Product Code	CSB-PA885732LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUV9
Immunogen	Recombinant Human GTPase IMAP family member 4 protein (5-101AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GTPase IMAP family member 4 (Immunity-associated nucleotide 1 protein) (IAN-1) (hIAN1) (Immunity-associated protein 4), GIMAP4, IAN1 IMAP4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	GIMAP4