



MDFIC Antibody, Biotin conjugated

Product Code	CSB-PA889174OD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P1T7
Immunogen	Peptide sequence from Human MyoD family inhibitor domain-containing protein (35-53AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MyoD family inhibitor domain-containing protein (I-mfa domain-containing protein) (hIC), MDFIC
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	MDFIC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.