



DICER1 Antibody, FITC conjugated

Product Code	CSB-PA891575LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UPY3
Immunogen	Recombinant Human Endoribonuclease Dicer protein (1721-1922AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoribonuclease Dicer (EC 3.1.26.3) (Helicase with RNase motif) (Helicase MOI), DICER1, DICER HERNA KIAA0928
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	DICER1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.