




MYH2 Antibody

Product Code	CSB-PA892458LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKX2
Immunogen	Recombinant Human Myosin-2 protein (1154-1300AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-2 (Myosin heavy chain 2) (Myosin heavy chain 2a) (MyHC-2a) (Myosin heavy chain IIa) (MyHC-IIa) (Myosin heavy chain, skeletal muscle, adult 2), MYH2, MYHSA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH2

Image



Western Blot

Positive WB detected in: HeLa whole cell lysate, Jurkat whole cell lysate

All lanes: MYH2 antibody at 1:1250

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 224, 80 kDa

Observed band size: 224 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.