



MYH2 Antibody, FITC conjugated

Product Code	CSB-PA892458LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKX2
Immunogen	Recombinant Human Myosin-2 protein (1154-1300AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Myosin-2 (Myosin heavy chain 2) (Myosin heavy chain 2a) (MyHC-2a) (Myosin heavy chain IIa) (MyHC-IIa) (Myosin heavy chain, skeletal muscle, adult 2), MYH2, MYHSA2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MYH2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.