



INPP5F Antibody, Biotin conjugated

Product Code	CSB-PA896692ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2H2
Immunogen	Peptide sequence from Human Phosphatidylinositide phosphatase SAC2 protein (187-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositide phosphatase SAC2 (EC 3.1.3.25) (Inositol polyphosphate 5-phosphatase F) (Sac domain-containing inositol phosphatase 2) (Sac domain-containing phosphoinositide 4-phosphatase 2) (hSAC2), INPP5F, KIAA0966 SAC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INPP5F
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.