



INPP5F Antibody, Biotin conjugated

Product Code	CSB-PA896692ND01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2H2
Immunogen	Peptide sequence from Human Phosphatidylinositide phosphatase SAC2 protein (187-204AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphatidylinositide phosphatase SAC2 (EC 3.1.3.25) (Inositol polyphosphate 5-phosphatase F) (Sac domain-containing inositol phosphatase 2) (Sac domain-containing phosphoinositide 4-phosphatase 2) (hSAC2), INPP5F, KIAA0966 SAC2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	INPP5F