



DYNC1LI1 Antibody, Biotin conjugated

Product Code	CSB-PA897312LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6G9
Immunogen	Recombinant Human Cytoplasmic dynein 1 light intermediate chain 1 protein (167-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytoplasmic dynein 1 light intermediate chain 1 (LIC1) (Dynein light chain A) (DLC-A) (Dynein light intermediate chain 1, cytosolic), DYNC1LI1, DNCLI1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	DYNC1LI1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.