



UTP18 Antibody, HRP conjugated

Product Code	CSB-PA897566LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5J1
Immunogen	Recombinant Human U3 small nucleolar RNA-associated protein 18 homolog protein (13-152AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	U3 small nucleolar RNA-associated protein 18 homolog (WD repeat-containing protein 50), UTP18, WDR50
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	UTP18
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.