




Recombinant *Macaca fascicularis* Interleukin-11 (IL11)

Product Code	CSB-EP324092MOV
Abbreviation	Recombinant Cynomolgus monkey IL11 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20808
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	<i>Macaca fascicularis</i> (Crab-eating macaque) (<i>Cynomolgus</i> monkey)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	PGPPPGSPRASDPRAELDSTVLLTRSLLEDTRQLTIQLKDKFPADGDHNLDS LPTLAMSAGALGALQLPSVLTRLRADLLSYLRHVQWLRRAMGSSLKTLEPELG TLQTRLDRLLRRLQLLMSRLALPQLPPDPPAPPLAPPSSTWGGIRAAHAILGGL HLTLDWAVRGLLLKTRL
Research Area	Cancer
Source	<i>E.coli</i>
Target Names	IL11
Protein Names	Recommended name: Interleukin-11 Short name= IL-11
Expression Region	22-199aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	25.4 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.