




Recombinant Escherichia coli O157:H7 Outer membrane protein X (ompX)

Product Code	CSB-EP359230EOD
Abbreviation	Recombinant E.coli O157:H7 ompX protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P0A919
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Escherichia coli O157:H7
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	ATSTVTGGY AQSDAQGQMNMKGGFNLKYRYEEDNSPLGVIGSFTYTEKSRTA SSGDYNKNQYYGITAGPAYRINDWASIYGVVGVGYGKFQTTEYPTYKHDTSD YGFSYGAGLQFNPMENVALDFSIEQSRIRSVDVGTWIAGVGYRF
Research Area	Others
Source	E.coli
Target Names	ompX
Protein Names	Recommended name: Outer membrane protein X
Expression Region	24-171aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	23.8 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.