




Recombinant Human papillomavirus type 16

Protein E7 (E7)

Product Code	CSB-YP365855HML
Abbreviation	Recombinant Human papillomavirus type 16 E7 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P03129
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Human papillomavirus type 16
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MHGDPTLTHEYMLDLQPETTDLYCYEQLNDSSEEEDEIDGPAGQAEPDRAHY NIVTFCKCDSTLRRCVQSTHVDIRLTLEDLLMGTLGIVCPICSQKP
Research Area	Epigenetics and Nuclear Signaling
Source	Yeast
Target Names	E7
Protein Names	Recommended name: Protein E7
Expression Region	1-98aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	13.0 kDa
Protein Length	Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.