




Recombinant Mouse Homeobox protein OTX2 (Otx2)

Product Code	CSB-EP017299MO
Abbreviation	Recombinant Mouse Otx2 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P80206
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MMSYLKQPPYAVNGLSLTTSMDLLHPSVGYPATPRKQRRERTTFTRAQLDV LEALFAKTRYPDIFMREEVALKINLPESRVQVWFKNRRRAKCRQQQQQQQNGG QNKVRPAKKKSSPAREVSSESGTSGQFTPPSSTSVPTIASSAPVSIWSPASIS PLSDPLSTSSSCMQRSYPMTYTQASGYSQGYAGSTSYFGGMDCGSYLTPMH HQLPGPGATLSPMGNTAVTSHLNQSPASLSTQGYGASSLGFNSTTDCLDYKD QTASWKLNFNADCLDYKDQTSSWKFQVL
Research Area	Developmental Biology
Source	E.coli
Target Names	Otx2
Protein Names	Recommended name: Homeobox protein OTX2 Alternative name(s): Orthodenticle homolog 2
Expression Region	1-289aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	37.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.