




Recombinant Candida albicans Hyphal wall protein 1 (HWP1), partial

Product Code	CSB-EP343212CZD
Abbreviation	Recombinant Candida albicans HWP1 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46593
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Candida albicans (strain SC5314 / ATCC MYA-2876) (Yeast)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GGGETEEALIQKRSYDYYQEPCDDYPQQQQQQEPCDYPQQQQQEEPCDYP QQQPQEPCDYPQQPQEPCDYPQQPQEPCDYPQQPQEPCDNPPQPDVPCD NPPQPDVPCDNPPQPDVPCDNPPQPDVPCDNPPQPDQPDNDNPPIPNIPTDWI PNIPTDWIPDIPEKPTTPATTPNIPA
Research Area	Neuroscience
Source	E.coli
Target Names	HWP1
Protein Names	Recommended name: Hyphal wall protein 1Alternative name(s): Cell elongation protein 2
Expression Region	27-203aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	24.2 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.