




Recombinant Mouse Cadherin-17 (Cdh17), partial

Product Code	CSB-EP882608MO
Abbreviation	Recombinant Mouse Cdh17 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9R100
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	INDVMYFQIDSKTGAISLTPEGSQELDPVKNPNSYNLVVSVKDMGGQSENSFSD TTYVDISIRENIWKAPEPVEIRENSTDP
Research Area	Signal Transduction
Source	E.coli
Target Names	Cdh17
Protein Names	Recommended name: Cadherin-17 Alternative name(s): BILL-cadherin Liver-intestine cadherin Short name= LI-cadherin P130
Expression Region	173-253aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	13.1 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients,



storage temperature and the stability of the protein itself.
Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.