




Recombinant Human Collagen alpha-2 (VI) chain (COL6A2), partial

Product Code	CSB-EP005752HU
Abbreviation	Recombinant Human COL6A2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P12110
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	ELSFVFLTDGVTGNDSLHESAHSRKQNVVPTVLALGSDVDMDVLTTLTSLGDRRAAVFHEKDYDSLAPGFFDFIR
Research Area	Signal Transduction
Source	E.coli
Target Names	COL6A2
Protein Names	Recommended name: Collagen alpha-2(VI) chain
Expression Region	941-1016aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	13.4 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.