




Recombinant Human Atrial natriuretic peptide receptor 2 (NPR2), partial

Product Code	CSB-EP016024HUa0
Abbreviation	Recombinant Human NPR2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P20594
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	RNLT LAVVLPEHNLSYAWAWPRVGP AVALAVEALGRALPVDLRFVSSELEGA CSEYLAPLSAVDLKLYHDPDLLLGP GCVYPAASVARFASHWRLPLLTAGAVAS GFSAKNDHYRTLVRTGPSAPKLGEFVVTLHG HFNWTARAALLYLDARTDDRP HYFTIEGVFEALQGSNLSVQH QVYAREPGGPEQATHFIRANGRIVYICGPLEML HEILLQAQRENLTNGDYVFFYLDVFGESLRAGPTRATGRP WQDNRTREQAQA LREAFQTVLVITYREPPNPEYQEFQNRLLIRAREDFGVELGPSLMNLIAGCFYD GILLYAEVLNETIQEGGTREDGLRIVEKMQGRRYHGV TGLVVMDKNNDRETDF VLWAMGDLDSGDFQPAAHYSGA EKQIWWTGRPIPWVKGAPPSDNPPCAF DL DDPSCDKTPLSTLAI
Research Area	Cardiovascular
Source	E.coli
Target Names	NPR2
Expression Region	23-458aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	54.5 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.