




Recombinant Mouse Merlin (Nf2), partial

Product Code	CSB-EP015741MO1
Abbreviation	Recombinant Mouse Nf2 protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P46662
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MAGAIASRMSFSSSLKRKQPKTFTVRIVTMDAEMEFNCEMKWKGKDLFDLVCR TLGLRETWFFGLQYTIKDTVAVLKMDDKVLVDHDVSKKEEPVTFHFLAKFY PENAEELVQEITQHLFFLQVKKQILDEKVYCPPEASVLLASYAVQAKYGDYDPSVHK RGFLAQEELLPKRVINLYQMTPEMWEERITAWYAEHRGRARDEAEMEYLKIAQ DLEMYGVNYFTIRNKKGTELLLGVDALGLHIYDPENRLTPKISFPWNEIRNISYS DKEFTIKPLDKKIDVFKFNSSKLRVNLKILQLCIGNHDLFMRRRKADSLEVQQ
Research Area	Cancer
Source	E.coli
Target Names	Nf2
Protein Names	Recommended name: Merlin Alternative name(s): Moesin-ezrin-radixin-like protein Neurofibromin-2 Schwannomin
Expression Region	1-320aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-Trx-tagged
Mol. Weight	55.9 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.