




# Recombinant Escherichia coli Universal stress protein D (uspD)

<b>Product Code</b>	CSB-EP364467ENV
<b>Abbreviation</b>	Recombinant E.coli uspD protein
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	P0AAB8
<b>Storage Buffer</b>	Tris-based buffer,50% glycerol
<b>Product Type</b>	Recombinant Proteins
<b>Immunogen Species</b>	Escherichia coli (strain K12)
<b>Purity</b>	Greater than 85% as determined by SDS-PAGE.
<b>Sequence</b>	MAYKHIGVAISGNEEDALLVNKALELARHND AHLTLIHIDDGLSELYPGIYFPAT EDILQLLKNKSDNKLYKLTKNIQWPCKLR IERGEMPETLLEIMQKEQC DLLVC GHHHSFINRLMPAYRGMINKMSADLLIVPFIDK
<b>Research Area</b>	Others
<b>Source</b>	E.coli
<b>Target Names</b>	uspD
<b>Protein Names</b>	Recommended name: Universal stress protein D
<b>Expression Region</b>	1-142aa
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	N-terminal 6xHis-tagged
<b>Mol. Weight</b>	20.4 kDa
<b>Protein Length</b>	Full Length

## Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.