



Recombinant Human Inhibitor of nuclear factor kappa-B kinase subunit beta (IKBKB)


Product Code	CSB-EP011572HU
Abbreviation	Recombinant Human IKBKB protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	O14920
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MSWSPSLTTQTCGAWEMKERLGTGGFGNVIRWHNQETGEQIAIKQCRQELS PRNRERWCLEIQIMRRLTHPNVVAARDVPEGMQNLAPNDLPLLAMEYCYGGD LRKYLNQFENCCGLREGAILTLLSDIASALRYLHENRIIHRDLKPENIVLQQGEQ RLIHKIIDLGYAKELDQGSLSCTSFVGTLLQYLAPELLEQQKYTVTVDYWSFGTLAF ECITGFRPFLPNWQPVQWHSKVRQKSEVDIVVSEDLNGTVKFSSSLPYPNNLN SVLAERLEKWLQLMLMWHPRQRGTDPTYGPNGCFKALDDILNLKLVHILNMVT GTIHTYPVTEDESLQSLKARIQQDTGIPEEDQELLQEAGLALIPDKPATQCISDG KLNIGHTLMDLVFLFDNSKITYETQISPRPQPESVSCILQEPKRNLAFFQLRK VWGQVWHSIQTLKEDCNRLQQGQRAAMNLLRNNSCLSKMKNMASMSQQ LKAKLDFFKTSIQIDLEKYSEQTEFGITSDKLLLAWREMEQAVELCGRENEVKLL VERMMALQTDIVDLQRSPMGRKQGGTLDLEEQARELYRRLREKPRDQRTE GDSQEMVRLLLQAIQSFEEKVRVIYTQLSKTVVCKQKALELLPKVEEVVSLMNE DEKTVVRLQEKRQKELWNLLKIACSKVRGPVSGSPDSMNASRLSQPGQLMSQ PSTASNSLPEPAKKSEELVAEAHNLCTLLENAIQDTVREQDQSFTALDWSWLQ TEEEHSCLEQAS
Research Area	Cancer
Source	E.coli
Target Names	IKBKB
Protein Names	Recommended name: Inhibitor of nuclear factor kappa-B kinase subunit beta Short name= I-kappa-B-kinase beta Short name= IKK-B Short name= IKK-beta Short name= IkbKB EC= 2.7.11.10 Alternative name(s): I-kappa-B ki
Expression Region	1-756aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged



Mol. Weight 92.5 kDa

Protein Length Full Length

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.