




Recombinant Mycobacterium tuberculosis Methionine aminopeptidase 2 (map)

Product Code	CSB-EP013716MVZ
Abbreviation	Recombinant Mycobacterium tuberculosis map protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P9WK18
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mycobacterium tuberculosis (strain CDC 1551 / Oshkosh)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MPSRTALSPGVLSPTRPVPNWIARPEYVGKPAAQEGSEPWWQTPEVIEKMRV AGRIAAGALAEAGKAVAPGVTTDELDRIAHEYLVDNGAYPSTLGYKGFPKSCC TSLNEVICHGIPDSTVITDGDIVNIDVTAYIGGVHGDNTATFPAGDVADEHRLLV DRTREATMRAINTVKPGRALSVIGRVIESYANRFGYNVVRDFTGHGIGTTFFHN GLVVLHYDQPAVETIMQPGMTFTIEPMINLGALDYEIWDDGWTVVTKDRKWA QFEHTLLVTDTGVEILTCL
Research Area	Signal Transduction
Source	E.coli
Target Names	map
Protein Names	Recommended name: Methionine aminopeptidase 2 Short name= MAP EC= 3.4.11.18Alternative name(s): Peptidase M
Expression Region	1-285aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged
Mol. Weight	36.9 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.