




Recombinant Mouse Apoptosis-associated speck-like protein containing a CARD (Pycard)

| | |
|--------------------------|---|
| Product Code | CSB-EP861664MO |
| Abbreviation | Recombinant Mouse Pycard protein |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q9EPB4 |
| Storage Buffer | Tris-based buffer,50% glycerol |
| Product Type | Recombinant Proteins |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | ≥ 85% as determined by SDS-PAGE. |
| Sequence | MGRARDAILDALENLSGDELKKFKMKLLTVQLREGYGRIPRGALLQMDAIDLTD KLVSYYLESYGLELTMVLRDMGLQELAEQLQTTKEESGAVAAAASVPAQSTA RTGHFVDQHRQALIARVTEVDGVLDALHGSVLTEGQYQAVRAETTSQDKMRK LFSFVPSWNLTKDSELLQALKEIHPYLVMDLEQS |
| Research Area | Cancer |
| Source | E.coli |
| Target Names | Pycard |
| Protein Names | Recommended name: Apoptosis-associated speck-like protein containing a CARD Short name= mASC Alternative name(s): PYD and CARD domain-containing protein |
| Expression Region | 1-193aa |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | N-terminal 10xHis-tagged and C-terminal Myc-tagged |
| Mol. Weight | 28.9 kDa |
| Protein Length | Full Length |
| Image | |



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.