




Recombinant Human Myc proto-oncogene protein (MYC), partial

Product Code	CSB-EP015270HU1
Abbreviation	Recombinant Human MYC protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P01106-1
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	SVCSTSSLYLQDL SAAASECIDPSVVFYPLNDSSSPKSCASQDSSAFSPSSD SLLSSTESSPQGSPEPLVLHEETPPTTSSDSEEEQEDEEEIDVVSVEKRQAPG KRSESGSPSAGGHSKPPHSPLVLKRCHVSTHQHNYYAAPSTRKDYPAAKRVK LDSVRVLRQISNNRKCTSPRSSDTEENVKRRTHNVLERQRRNELKRSFFALRD QIPELENNEKAPKVILKKATAYILSVQAEEQKLISEEDLLRKRREQLKHKLEQL RNSCA
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	MYC
Expression Region	169-439aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-GST-tagged
Mol. Weight	61.7 kDa
Protein Length	Partial of Isoform 1

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.