




Recombinant Bovine Cathelicidin-6 (CATHL6)

Product Code	CSB-YP345751BO
Abbreviation	Recombinant Bovine CATHL6 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P54228
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Bos taurus (Bovine)
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	GRFKRFRKKFKKLFKKLSPVIPLLHLG
Research Area	Others
Source	Yeast
Target Names	CATHL6
Protein Names	Recommended name: Cathelicidin-6 Alternative name(s): Antibacterial peptide BMAP-27 Myeloid antibacterial peptide 27
Expression Region	132-158aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag-Free
Mol. Weight	3.3 kDa
Protein Length	Full Length of Mature Protein

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.