




Recombinant Rat COMM domain-containing protein 5 (Commd5)

Product Code	CSB-EP880645RA
Abbreviation	Recombinant Rat Commd5 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9ERR2
Storage Buffer	Tris-based buffer,50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	SALGAAAPYLHHPADSHSGRVSFLGSQPSPEVTAVAQLLKDLDRSTFRKLLKL VVGALHGKDCREAVEQLGASANLSEERLAVLLAGTHTLLQQALRLPPASLKPD AFQEELQELGIPQDLIGDLASLAFGSQRPLLDSPAQQQGSSLPHVSYFRWRVD VAISTSAQSRSLQPSVLMQLKLTGSAHRFEVPIAKFQELRYSVALVLKEMAEL EKKCERKLQD
Research Area	Epigenetics and Nuclear Signaling
Source	E.coli
Target Names	Commd5
Protein Names	Recommended name: COMM domain-containing protein 5Alternative name(s): Hypertension-related calcium-regulated gene protein Short name= HCaRG
Expression Region	2-224aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	31.8 kDa
Protein Length	Full Length of Mature Protein
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.