




Recombinant Streptomyces viridosporus Probable hydrogen peroxide-inducible genes activator (oxyR)

Product Code	CSB-EP895843SNAG
Abbreviation	Recombinant Streptomyces viridosporus oxyR protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q9X5P2
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Streptomyces viridosporus
Purity	≥ 90% as determined by SDS-PAGE.
Sequence	MRKRRRQPSLAQLRAFAAVA EHLHFRDAAAAIGMSQPALSGAVSALEESLGV TLLERTTRKVLLSPAGERLAARAKSVLA EVGALVEEADALQAPFTGVLRLGVIP TVAPYVLPTVLR LVH DRYPRLDLQVHEEQ TASLLDGLTGGRLDLLLAVPLGVP GVVELPLFDEDFVLVTPLEHGLGGREGIPRKALRELNLLLLDEGHCLRDQALDI CREAGSAGVAATTTAAGLSTLVQLVAGGLGV TLLPHTAVQVETTRSGRLLTGR FADPAPGRRIALAMRTGAARAAEYEELAAALREAMRDL PVRIVRD
Research Area	Others
Source	E.coli
Target Names	oxyR
Protein Names	Recommended name: Probable hydrogen peroxide-inducible genes activator
Expression Region	1-312aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	40.6 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.