




Recombinant Mouse Receptor-interacting serine/threonine-protein kinase 1 (Ripk1)

Product Code	CSB-EP720181MOa2
Abbreviation	Recombinant Mouse Ripk1 protein
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q60855
Storage Buffer	Tris-based buffer, 50% glycerol
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	MQPDMSLDNIKMASSDLLEKTDLDSSGGFGKVSCLCYHRSHGFVILKKVYTGPNR AEYNEVLLEEGKMMHRLRHSRVVKLLGIIIEEGNYSVMMEYMEKGNLMHVLKT QIDVPLSLKGRIIVEAIEGMCYLHDKGVIHKDLKPENILVDRDFHIKIADLGVASFK TWSKLTKEKDNKQKEVSSTTKNNGGTLYYMAPEHLNDINAKPTEKSDVYSF GIVLWAIFAKKEPYENVICTEQFVICIKSGNRPNVEEILEYCPREIISLMERCWQA IPEDRPTFLGIEEEFRPFYLSHFEEYVEEDVASLKKEYPDQSPVLQRMFSLQHD CVPLPPSRNSEQPGSLHSSQGLQMGPVEESWFSSSPEYPQDENDRSVQAK LQEEASYHAFGIFAQKQTKPQPRQNEAYNREEERKRRVSHDPFAQQRARENI KSAGARGHSDPSTTSRGIAVQQLSWPATQTVWNNGLYNQHGFGTGTGVVWY PPNLSQMYSTYKTPVPETNIPGSTPTMPYFSGPVADDLIKYTIFNSSGIQIGNHN YMDVGLNSQPPNNTCKEESTSRHQAIFDNTTSLTDEHLNPIRENLGRQWKNC ARKLGFTESQIDEIDHDYERDGLKEKVYQMLQKWLMMREGTKGATVGKLAQAL HQCCRIDLLNHLIRASQS
Research Area	Cancer
Source	E.coli
Target Names	Ripk1
Expression Region	1-656aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-SUMO-tagged
Mol. Weight	87.8 kDa
Protein Length	Full Length
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.