




Recombinant Mouse Phospholipase A-2-activating protein (Plaa), partial

Product Code	CSB-BP018107MO
Abbreviation	Recombinant Mouse Plaa protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P27612
Product Type	Recombinant Proteins
Immunogen Species	Mus musculus (Mouse)
Purity	≥ 85% as determined by SDS-PAGE.
Sequence	TGAGRYMPGSAGMDTTMTGVDPFTGNSAYRSAASKTVNIYFPKKEALTFDQA NPTQILGKLKELNGTAPEEKKLTEDDLVLLEKILSLIC
Research Area	Immunology
Source	Baculovirus
Target Names	Plaa
Protein Names	Recommended name: Phospholipase A-2-activating protein Short name= PLA2P Short name= PLAP
Expression Region	495-584aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 10xHis-tagged and C-terminal Myc-tagged
Mol. Weight	13.6 kDa
Protein Length	Partial

Image



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.