




Recombinant Human Follicle-stimulating hormone receptor (FSHR), partial

Product Code	CSB-MP009021HU
Abbreviation	Recombinant Human FSHR protein, partial
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P23945
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	Greater than 85% as determined by SDS-PAGE.
Sequence	CHHRICHCSNRVFLCQESKVTEIPSDLPRNAIELRFVLTCLRVIQKGAFGSGFDL EKIEISQNDVLEVIEADVFSNLPKLHEIRIEKANNLLYINPEAFQNLPNLQYLLISN TGIKHLPDVHKKIHSLSLQKVLDDIQDNINIHTIERNSEFVGLSFESVILWLNKNGIQEIH NCAFNGTQLDELNLSNNDNLEELPNDVFHGGASGPVILDISRTRIHSLPSYGLN LKKLRARSTYNLKKLPTLEKLVALMEASLTYPHCCAFANWRRQISELHPICNK SILRQEVDMYMTQARGQRSSLAEDNESSYRSGFDMTYTEFDYDLCNEVVDVTC SPKPDFAFNPCEDIMGYNILR
Research Area	Neuroscience
Source	Mammalian cell
Target Names	FSHR
Protein Names	Recommended name: Follicle-stimulating hormone receptor Short name= FSH- R Alternative name(s): Follitropin receptor
Expression Region	18-366aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	N-terminal 6xHis-tagged
Mol. Weight	43.5 kDa
Protein Length	Partial
Image	



(Tris-Glycine gel) Discontinuous SDS-PAGE (reduced) with 5% enrichment gel and 15% separation gel.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.