



Ace Antibody, HRP conjugated

Product Code	CSB-PA001149LB01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P47820
Immunogen	Recombinant Rat Angiotensin-converting enzyme protein (1046-1263AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Angiotensin-converting enzyme (ACE) (EC 3.2.1.-) (EC 3.4.15.1) (Dipeptidyl carboxypeptidase I) (Kininase II) (CD antigen CD143) [Cleaved into: Angiotensin-converting enzyme, soluble form], Ace, Dcp1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Ace
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.