



Ace Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA001149LB01RA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | P47820 |
| Immunogen | Recombinant Rat Angiotensin-converting enzyme protein (1046-1263AA) |
| Raised In | Rabbit |
| Species Reactivity | Rat |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Angiotensin-converting enzyme (ACE) (EC 3.2.1.-) (EC 3.4.15.1) (Dipeptidyl carboxypeptidase I) (Kininase II) (CD antigen CD143) [Cleaved into: Angiotensin-converting enzyme, soluble form], Ace, Dcp1 |
| Immunogen Species | Rattus norvegicus (Rat) |
| Research Area | Cardiovascular |
| Target Names | Ace |