



# Ace Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA001149LD01RA
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P47820
<b>Immunogen</b>	Recombinant Rat Angiotensin-converting enzyme protein (1046-1263AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Rat
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Angiotensin-converting enzyme (ACE) (EC 3.2.1.-) (EC 3.4.15.1) (Dipeptidyl carboxypeptidase I) (Kininase II) (CD antigen CD143) [Cleaved into: Angiotensin-converting enzyme, soluble form], Ace, Dcp1
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Research Area</b>	Cardiovascular
<b>Target Names</b>	Ace
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.