



CEACAM5 Antibody, Biotin conjugated

Product Code	CSB-PA005165LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P06731
Immunogen	Recombinant Human Carcinoembryonic antigen-related cell adhesion molecule 5 protein (1-255AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carcinoembryonic antigen-related cell adhesion molecule 5 (Carcinoembryonic antigen) (CEA) (Meconium antigen 100) (CD antigen CD66e), CEACAM5, CEA
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CEACAM5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.