



Gsta1 Antibody, Biotin conjugated

Product Code	CSB-PA009970HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13745
Immunogen	Recombinant Mouse Glutathione S-transferase A1 protein (2-223AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Glutathione S-transferase A1 (EC 2.5.1.18) (GST class-alpha member 1) (Glutathione S-transferase Ya) (Glutathione S-transferase Ya1) [Cleaved into: Glutathione S-transferase A1, N-terminally processed], Gsta1, Gsta Gsty
Immunogen Species	Mus musculus (Mouse)
Research Area	Signal Transduction
Target Names	Gsta1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.