



HAP1 Antibody, HRP conjugated

Product Code	CSB-PA010129LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54257
Immunogen	Recombinant Human Huntingtin-associated protein 1 protein (328-580AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Huntingtin-associated protein 1 (HAP-1) (Neuroan 1), HAP1, HAP2 HLP1
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	HAP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.