



NEURL1B Antibody, Biotin conjugated

Product Code	CSB-PA015722LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8MQ27
Immunogen	Recombinant Human E3 ubiquitin-protein ligase NEURL1B protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase NEURL1B (EC 2.3.2.27) (Neutralized-2) (NEUR2) (Neutralized-like protein 1B) (Neutralized-like protein 3) (RING-type E3 ubiquitin transferase NEURL1B), NEURL1B, NEURL3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NEURL1B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.