



# NEURL1B Antibody, Biotin conjugated

Product Code	CSB-PA015722LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	A8MQ27
Immunogen	Recombinant Human E3 ubiquitin-protein ligase NEURL1B protein (1-275AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase NEURL1B (EC 2.3.2.27) (Neuralized-2) (NEUR2) (Neuralized-like protein 1B) (Neuralized-like protein 3) (RING-type E3 ubiquitin transferase NEURL1B), NEURL1B, NEURL3
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NEURL1B